

PRODUCT INFORMATION

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806 1766 Clone ID

GFAP ALXDRD Synonyme Host Species Rabbit

Biotinylated Anti-GFAP(68-377) antibody(17E6), Rabbit mAb Description

Delivery 2-3 weeks Uniprot ID P14136 lgG type Rabbit IgG Clonality Monoclonal Reactivity Human Application ELISA

Recommended Dilutions ELISA 1:5000-10000

Purified from cell culture supernatant by affinity chromatography Purification

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specific instructions or reconstitution. Store at 2-90° C 0-80°°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct isoforms, [provided by RefSeq. 0ct 2008]

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only Conjugate

Notice that the complete of th DIMA Disclaimer

