Cat. No. DMC100478B



PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806 bMC478 Clone ID

GPR87

FKSG78; GPR95; KPG_002

Host Species Rabbit

Biotinylated Anti-GPR87 antibody(DMC478); IgG1 Chimeric mAb Description

Uniprot ID O9BY21

IgG type Rabbit/Human Fc chimeric IgG1

Reactivity Human Application Flow Cyt Recomme Dilutions Flow Cyt 1:100 Delivery

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-couple receptor genes on chronosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589). [provided by RefSeq: Nov 2011] Storage & Shipping

Background

Usage Research use only

Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer




Figure 1. Flow cytometry analysis with $15\mu g/mL$ Biotinylated Anti-GPR87 (DMC478) on HEK293 cells transfected with human GPR87 (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

