

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\https://www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806
HMBDD01:HMBD-001:HMBD 001:HMB 001 Common Name

Conjugate

HER3;FERLK;LCCS2;VSCN1;ErbB-3;c-erbB3;erbB3-S;MDA-BF-1;c-erbB-3;p180-ErbB3;p45-sErbB3;p85-sErbB3 Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species lgG1 IgG type Reactivity Target HER3 Uniprot ID

Biotinylated Anti-Her3(nezutatug biosimilar) mAb Description

Delivery

6-8 weeks
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

