Cat. No. DME101031B



## **PRODUCT INFORMATION**

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products/public/class-woocommerce-print-products-public.php on line 2806
211. Clone ID

Synonyme Host Species Rabbit

Biotinylated Anti-Human IgG antibody(2H2), Rabbit mAb Description

Delivery 2-3 weeks Uniprot ID N/A IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Applications ELISA Recommended Dilutions ELISA 1:5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Access Ac Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Usage Research use only

Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

> Email: info@dimabio.com Website: www.dimabio.com

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

