Cat. No. DMC100367B



PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in Cs\www.root\wirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806 bMC367 Clone ID

ICAM1;BB2;CD54;P3.58 Synonyme

Host Species Rabbit

Biotinylated Anti-ICAM-1 antibody(DMC367); IgG1 Chimeric mAb Description

Delivery 2-3 weeks P05362 Uniprot ID

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Applications Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Leguencing or reverse specinic instructions of reconstitution.

Store at 20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freeign and thawing). Lyophilized proteins are shipped at ambient temperature.

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD1a: CD18; or CD11b: CD18 and is also exploited by Rhinovirus as a receptor.

Research use only Storage & Shipping Background

> Email: info@dimabio.com Website: www.dimabio.com

Conjugate Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

