Cat. No. DMC100422B



## PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806 bMC422 Clone ID

IBP-7; TAF; IGFBP-rP1; MAC25 protein Synonyme

Host Species Rabbit

Biotinylated Anti-IGFBP7 antibody(DMC422); IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID Q16270

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Application Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

specific instructions of reconstitution.

Store at 20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature, repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature, repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature, and the store of the st

Usage Coniugate

Background

protein s. All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer



Email: info@dimabio.com Website: www.dimabio.com

