

PRODUCT INFORMATION

Clone ID

II 18RP FVH;IL18BPa Synonyme Host Species Rabbit

Biotinylated Anti-IL18BP antibody(1G1), IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID 095998

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Application ELISA

Recommended Dilutions ELISA 1/5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specific instructions of reconstitution.

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

The protein encoded by this gene functions as an inhibitor of the promitianmentory of the left of the protein for the protein f Background

Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate Research use only Biotinylated

DIMA Disclaimer

.ein set, All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

