Cat. No. DME100218B



PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bM218 Clone ID

II 1884

CD218a; CDw218a; IL-1Rrp; IL1RRP; IL18RA Synonyme

Host Species Rabbit

Biotinylated Anti-IL18RA antibody(DM218); Rabbit mAb Description Delivery 2-3 weeks

Q13478 Uniprot ID lgG type Rabbit IgG Clonality Monoclonal Reactivity Human Application ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Specific instructions or reconstitution. Store at 20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at .80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

The protein recorded by this gene is a cytokine receptor that belongs to the interference control of the protein control of the protein recorded by this gene is a cytokine receptor that belongs to the interference control of the protein specifically binds interference in the protein specifically binds in the p Background

> Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate

, roten s All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

