

## **PRODUCT INFORMATION**

Warning. Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806 bMC271. Clone ID

IL-1; IL1-BETA; IL1F2 Synonyme

Host Species Rabbit

Biotinylated Anti-IL1B antibody(DMC271); IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID P01584

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Applications Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specific instructions of reconstitution.

Store at .20°C to .80°C for 12 months in Jophilized form. After reconstitution in the interest of analysis for specific instructions of reconstitution in Jophilized form. After reconstitution in the interest of t Background

> Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate

Research use only

Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. **DIMA Disclaimer** 

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

