

PRODUCT INFORMATION

Warning. Undefined variable ShasAttributeValueDescription in Chwwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806
[IM 33.L(M33.)] Common Name

Conjugate

Biotinylated
CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384; zcytoR17 Synonyms

Applications ELISA, Flow Cyt

Recommended Dilutions ELISA 1:5000-10000, Flow Cyt 1:100

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species

IgG2(C14S,R16K,E20G,S21G,C102S,H147Q,R234Q,S257A,Q298E) IgG type

Reactivity Target IL31RA Uniprot ID

Biotinylated Anti-IL31RA(nemolizumab biosimilar) mAb Description

Delivery

Jet (Areid 2-3 weeks
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. Storage & Shipping

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

