Cat. No. DMC100278B



PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in Cs\www.root\wirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.chass-woocommerce-print-products-public.php on line 2806
DMC278 Clone ID

JAM-A CD321; JAM;JAM1; JAMA; JCAM; KAT; PAM-1 Synonyme

Host Species Rabbit

Biotinylated Anti-JAM-A antibody(DMC278); IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID Q9Y624

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Application Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

specific instructions of reconstitution.

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid respected freezing and chaving). Lyophilized proteins are shipped at ambient temperature.

For example, the protein store of the protein store shipped at ambient temperature.

For example, the protein continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. The protein encoded by this immunoglobulin superfamily gene member is an important regulator of titing function assembly in epithelia. In addition, the encoded protein can act as (1) a receptor for recovings; (2) a ligand for the integrin LFA1; involved in leukocyte transmigration; and (3) a platelet receptor. Multiple 5' alternatively spliced variants; encoding the same protein; have been identified but their biological validity has not been established.

Research use only

Usage Coniugate

Background

, rotein se All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

> Email: info@dimabio.com Website: www.dimabio.com

