

PRODUCT INFORMATION

Warning: Undefined variable ShacktributeValueDescription in CAlwawroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public.php on line 2806

Conjugate

kidney associated antigen 1;kidney-associated antigen 1;MGC78738;RU2 antisense gene protein;RU2AS Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species lgG1 IgG type Reactivity Target KAAG1 Uniprot ID

Biotinylated Anti-KAAG1(olintatug biosimilar) mAb Description

Delivery

6.8 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage



