

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 IPH 41;IPH 41;IPH 41;IPH 41;IPH 41;IPH 41;IPH 41;IPH 41;IPH 41;IPH 41
Conjugate	Biotinylated
Synonyms	3DL2;p140;NKAT4;CD158K;NKAT-4;NKAT4B;KIR-3DL2
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	KIR3DL2
Uniprot ID	P43630
Description	Biotinylated Anti-KIR3DL2(lacutamab biosimilar) mAb
Delivery	6-8 weeks
Yefei Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

