

## **PRODUCT INFORMATION**

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806
bMI84 Clone ID

CLEC15A; MAFA; MAFAL Synonyme

Host Species Rabbit

Biotinylated Anti-KLRG1 antibody(DM184); Rabbit mAb Description

Delivery 2-3 weeks Uniprot ID Q96E93 IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Application ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Store at -20°C to -80°C for 12° months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Natural killer (MK) cells are lymphocytes that can mediate lysis of certain tumor cells and virune-infected cells without previous activation. They can also Natural killer (MK) cells are lymphocytes that can mediate lysis of certain tumor cells and virune-infected cells without previous activation. They can also have been considered to the control of the cells of the Background

Usage Conjugate Research use only Biotinylated

, rotein se All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

