

## **PRODUCT INFORMATION**

Warning: Undefined variable ShasAttributeValueDescription in C\u00e4wwroot\u00fcmirror.dimabio.com\u00e4wp-content\plugins\u00e4woocommerce-print-products-public.php on line 2806 3811 Clone ID

MUC16 Synonyme Host Species Rabbit

Biotinylated Anti-MUC16 antibody(3B11), IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID Q8WXI7

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Application Flow Cyt Recommended Dilutions Flow Cyt 1/100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 26% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

specific instructions of reconstitution.

Store at .20°C to .80°C for 12 months in lyophilized form. After reconstitution, in not intended for use within a month, aliquot and store at -98°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

The interpolation of the control o Background

Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate

Research use only

Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. **DIMA Disclaimer** 



