

## **PRODUCT INFORMATION**

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-bold
618 Clone ID

CMT1F; CMT2E; CMTDIG; NF-L; NF68; NFL; PPP1R110 Synonyme

Host Species Rabbit

Biotinylated Anti-NEFL(89-400) antibody(6D8), Rabbit mAb Description

Delivery 2-3 weeks P07196 Uniprot ID lgG type Rabbit IgG Clonality Monoclonal Reactivity Human Applicatio ELISA

Recommende Dilutions ELISA 1:5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specific instructions or reconstructions. Some at 20°C to 80°C for 12°C to 80°C for 12°C to 80°C for 12°C months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Neurofilaments are type (I) intermediate filament heteropolymers composed of light, medium, and heavy chains. Neurofilaments comprise the Neurofilaments are typed in the state of Background

Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate Research use only Biotinylated

de seç All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

