

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 AEX-1188;AEX1188;AEX 1188;INCMGA-0012;INCMGA 0012;INCMGA0012;INCMGA12;INCMGA-12;INCMGA 12;MGA-012;MGA012;MGA 012;MGA-12;MGA12;MGA 12;Zynyz
Conjugate	Biotinylated
Synonyms	PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG4
Reactivity	Human
Target	PD-1
Uniprot ID	Q15116
Description	Biotinylated Anti-PD-1(retifanlimab biosimilar) mAb
Delivery	6-8 weeks
Yefei_Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

