Cat. No. BME100703B



## **PRODUCT INFORMATION**

Warning: Undaffixed variable ShasAttributeValueDescription in CAywwroottmirror.dimable.com/wp-content/plugins/woocommerce-picks-gov/thecay-pitht-class-woocommerce-pitht-products-public-pith on line 2006.
CHS-007/CHS-007/CHS-007/CHS-007/CHS-0015/SD-015/SD **Common Name** 

Conjugate

Synonyms Applications PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1

Recommended Dilutions N/A

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Host Species** 

N/A Reactivity Human Target PD-1 Uniprot ID Q15116

Biotinylated Anti-PD-1(toripalimab biosimilar) mAb

Delivery 6-8 weeks

To weeks. Store at 20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature and thawing at the store at -80°C (Avoid Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background

Usage Research use only

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

