

## PRODUCT INFORMATION

<b>Common Name</b>	<b>Warning:</b> Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 STI-1014;STI1014;STI 1014;STI-A1014;STI A1014;ZKAB-001;ZKAB001;ZKAB 001;ZKAB-1;ZKAB1;ZKAB 1
<b>Conjugate</b>	Biotinylated
<b>Synonyms</b>	B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	N/A
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	PD-L1
<b>Uniprot ID</b>	Q9NZQ7
<b>Description</b>	Biotinylated Anti-PD-L1(socazolimab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Yefei Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

