Cat. No. DMC100428B



PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in Cs\www.root\wirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.chass-woocommerce-print-products-public.php on line 2806
bMC428 Clone ID

SSTR2 SRIF-1; SS2R Synonyme Host Species Rabbit

Biotinylated Anti-SSTR2 antibody(DMC428); IgG1 Chimeric mAb Description

Delivery 2-3 weeks P30874 Uniprot ID

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Applications Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specinic instructions or reconstitutions. Store at 2-02°C to -80°C for 12°C to -80°C for 12°C for 12°C

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

Conjugate

Notice that the complete of th Julian se, DIMA Disclaimer

