

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 DSTP-3086S, DSTP3086S, MSTP2109A, RG-7450, RG7450
Conjugate	Biotinylated
Synonyms	PRSS24, STEAP
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 15µg/ml
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	STEAP1
Uniprot ID	Q9UHE8
Description	Biotinylated Anti-STEAP1(Vandortuzumab biosimilar) mAb
Delivery	In Stock
Yefei Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 15 µg/ml Biotinylated Anti-STEAP1(Vandortuzumab biosimilar) mAb on SNU-5 cell line (Blue histogram) or 293T (Red histogram).

