

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bM178 Clone ID

Synonyme

Host Species Rabbit

Biotinylated Anti-TIGIT antibody(DM178); Rabbit mAb Description

Delivery 2-3 weeks Uniprot ID Q495A1 IgG type Rabbit IgG Clonality Monoclonal Reactivity Applications ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specinic instructions or reconstitutions. Store at 2-70°C to -80°C for 12°C to -80°C for 12°C for 12°C

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only Conjugate

Nescentro Use Unity
Biotinylated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. Azin se. DIMA Disclaimer

