

## **PRODUCT INFORMATION**

Warning: Undefined variable ShacktributeValueDescription in CAwwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\upballc\class-woocommerce-print-products-public.php on line 2806 TSR-022 Common Name

Conjugate

Biotinylated HAVCR2;TIM3;TIMD3;FLJ14428;KIM3 Synonyms

Applications ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species IgG type IgG4

Reactivity Human Target TIM3 Uniprot ID Q8TDQ0

Biotinylated Anti-TIM3 (cobolimab biosimilar) mAb Description

Delivery

2-3 weeks
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Background Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

