

## **PRODUCT INFORMATION**

Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 AGX.101 Common Name

Biotinylated M3S1;TAAL6 Conjugate Synonyms Applications Flow Cyt Recommended Dilutions Flow Cyt 1:100

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species

IgG type lgG1 Reactivity Target TM4SF1 Uniprot ID

Biotinylated Anti-TM4SF1 (AGX101 biosimilar) mAb Description

Delivery

2-3 weeks
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

