

PRODUCT INFORMATION

Common Name	GEN1029;GEN-1029;Hx-DR5-01;Hx-DR5-05;iDD 004;iDD-004;iDD004
Synonyms	DR5;CD262;KILLER;TRICK2;TRICKB;ZTNFR9;TRAILR2;TRICK2A;TRICK2B;TRAIL-R2;KILLER/DR5
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	TNFRSF10B
Uniprot ID	O14763
Description	Biotinylated Anti-TNFRSF10B(benufutamab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Biotinylated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

DIMABIO CONFIDENTIAL

