

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C\u00e4wwroot\u00fcmirror.dimabio.com\u00e4wp-content\plugins\u00e4woocommerce-print-products-public.php on line 2806 AMG162 Common Name

Conjugate

AMG102
Biotinylated
CD254;hRANKL2;ODF;OPGL;OPTB2;RANKL;sOdf;TNLG6B;TRANCE Synonyms

Applications Flow Cyt Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species IgG type IgG2 Reactivity Target TNFSF11 Uniprot ID

Biotinylated Anti-TNFSF11 (denosumab biosimilar)mAb Description

Delivery

. *C (wost 2-3 weeks
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

