

## **PRODUCT INFORMATION**

Warning: Undefined variable ShasAttributeValueDescription in C\\u00edword\u00edmirror.dimabio.com\u00edwp-content\plugins\u00edwoocommerce-print-products\u00edpublic\u00edlass-woocommerce-print-products-public.php on line 2806
AMG 157, AMG-157, MCD-1929, MCD19292, anti-TSP-Receptor ITSP-R: Cri2D, tezepelumab-ekko Common Name

Conjugate

Biotinylated
Thymic stromal lymphopoietin Synonyms

Applications ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species

IgG type IgG2 Reactivity Human Target TSLP Uniprot ID

Biotinylated Anti-TSLP(tezepelumab biosimilar) mAb Description

Delivery

C local 2-3 weeks
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

