

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\unwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806
244 Clone ID

LICHI 1

HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1 Synonyme

Host Species

Biotinylated Anti-UCHL1 antibody(2A4), Rabbit mAb Description

Delivery 2-3 weeks P09936 Uniprot ID IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Application ELISA; WB

Recommended Dilutions ELISA 1:5000-10000; WB 1:1000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Store at 20°C to 80°C for 12° months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol protease that hydrolyzes a peptide bond at the C-terminal givcine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease [provided by RefSeq., Sep 2009]

Usage Research use only Conjugate

Notice that the complete of the control of the cont ein set, DIMA Disclaimer

> Email: info@dimabio.com Website: www.dimabio.com

