

PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in Cs\www.root\wirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806
bMC389 Clone ID

Synonyme Host Species Rabbit

Biotinylated Anti-VWF antibody(DMC389); IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID P04275

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Applications Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Store at 20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. This gene recodes a glycoprotein involved in hemostasis. The encoded preproportien is proteolytically processed following assembly into large that the processes of the control of the

Email: info@dimabio.com Website: www.dimabio.com

Background

Usage Conjugate Research use only Biotinylated

æin seg All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

