

PRODUCT INFORMATION

PD-1 **Target**

PD1; PDCD1; CD279; SLEB2; hPD-1; hPD-l; hSLE1 **Synonyms**

Recombinant Canine PD-1 protein with C-terminal **Description**

10×His tag

Delivery In Stock

Uniprot ID NP_001301026.1

Expression Host HEK293

Tag C-10×His tag

Molecular

Background

PD-1(Leu25-Gly168) 10×His tag Characterization

The protein has a predicted molecular mass of **Molecular Weight**

17.5 kDa after removal of the signal peptide. The apparent molecular mass of dPD-1-His is

approximately 15-35 kDa due to glycosylation.

The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue

Purity

staining.

Lyophilized from sterile PBS, pH 7.4. Normally 5 % Formulation & Reconstitution

- 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

temperature.

Programmed cell death protein 1 (PDCD1) is an

immune-inhibitory receptor expressed in activated T cells; it is involved in the regulation of T-cell functions, including those of effector CD8+

T cells. In addition, this protein can also promote the differentiation of CD4+ T cells into T regulatory cells. PDCD1 is expressed in many

types of tumors including melanomas, and has demonstrated to play a role in anti-tumor immunity. Moreover, this protein has been shown

to be involved in safeguarding against

autoimmunity, however, it can also contribute to the inhibition of effective anti-tumor and antimicrobial immunity. [provided by RefSeq, Aug

> Email: info@dimabio.com Website: www.dimabio.com

20201

Usage Research use only

Conjugate Unconjugated








Figure 1. Canine PD-1 Protein, His Tag on SDS-PAGE under reducing condition.

Canine PD-1, His Tagged protein ELISA

0.2 μg of Canine PD-1, His tagged protein per well




Figure 2. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Canine PD-1 Protein, His Tag (PME-D100005) can bind Anti-Canine PD-1(INTERVET 4F12) mAb (BME100202) in a linear range of 0.64–80 ng/mL.

Email: info@dimabio.com Website: www.dimabio.com



Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)






Figure 3. Flow cytometry analysis with $15\mu g/mL$ Canine PD-1 Protein, His Tag(PME-D100005) on HEK293 cells transfected with Canine PD-L1 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Email: info@dimabio.com Website: www.dimabio.com

