

PRODUCT INFORMATION

Target	EGFR
Description	Monoclonal Cell Line Derived from K562 Cells, Engineered for Stable Expression of Human EGFR Using Lentiviral Technology
Host Cells	K562
Uniprot ID	P00533
Applications	FACS Data
Growth media	RPMI-1640+10% FBS+1% P.S+2 ug/mL Puromycin
Package	5E6 Cells/mL
Host Species	Human
Suggested Control	SKU: BME100034
Warranty and Disclaimer	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
Storage & Shipping	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
Synonyms	EGFR;ERBB;ERBB1;HER1;PIG61;mENA
Background	The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer.
Usage	For research use only.



Hu_EGFR K562 Cell Line




Figure 1. Flow cytometry analysis of human EGFR overexpression using Hu_EGFR K562 Cell Line (Cat. No. CEL100005) and Anti-EGFR (Cetuximab biosimilar) mAb (Cat. No. BME100034)

