


**PRODUCT INFORMATION**

<b>Target</b>	CD112
<b>Description</b>	Monoclonal Cell Line Derived from CHO-S Cells, Engineered for Stable Expression of Human CD112 Using Lentiviral Technology
<b>Host Cells</b>	CHO-S
<b>Uniprot ID</b>	Q92692
<b>Applications</b>	FACS Data
<b>Growth media</b>	DMEM+10% FBS+1% P.S+Gln+2 ug/mL Puromycin
<b>Package</b>	5E6 Cells/mL
<b>Host Species</b>	Human
<b>Suggested Control</b>	SKU: DMC100224
<b>Warranty and Disclaimer</b>	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
<b>Storage &amp; Shipping</b>	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
<b>Synonyms</b>	NECTIN2; HVEB; PRR2; PVRL2; PVRR2
<b>Background</b>	This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus; and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants; encoding different isoforms; have been characterized.
<b>Usage</b>	For research use only.



## Hu\_CD112 CHO-S Cell Line



- Human IgG
- Anti-CD112 antibody(DMC224);  
IgG1 Chimeric mAb (SKU: DMC100224)

Figure 1. Flow cytometry analysis of human CD112 overexpression using Hu\_CD112 CHO-S Cell Line (Cat. No. CEL100077) and Anti-CD112 antibody(DMC224)IgG1 Chimeric mAb (Cat. No. DMC100224)

DIMABIO CONFIDENTIAL

