


PRODUCT INFORMATION

Target	CD79B
Description	Monoclonal Cell Line Derived from CHO-S Cells, Engineered for Stable Expression of Human CD79B Using Lentiviral Technology
Host Cells	CHO-S
Uniprot ID	P40259
Applications	FACS Data
Growth media	DMEM+10% FBS+1% P.S+Gln+2 ug/mL Puromycin
Package	5E6 Cells/mL
Host Species	Human
Suggested Control	SKU: BME100171
Warranty and Disclaimer	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
Yefei_Storage	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
Synonyms	AGM6;B29;IGB
Background	The B lymphocyte antigen receptor is a multimeric complex that includes the antigen-specific component, surface immunoglobulin (Ig). Surface Ig non-covalently associates with two other proteins, Ig-alpha and Ig-beta, which are necessary for expression and function of the B-cell antigen receptor. This gene encodes the Ig-beta protein of the B-cell antigen component. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]
Usage	For research use only.



Hu_CD79B CHO-S Cell Line



-  Human IgG
-  Anti-CD79B(polatuzumab biosimilar) mAb (SKU: BME100171)

Figure 1. Flow cytometry analysis of human CD79B overexpression using Hu_CD79B CHO-S Cell Line (Cat. No. CEL100079) and Anti-CD79B(polatuzumab biosimilar) mAb (Cat. No. BME100171)

DIMABIO CONFIDENTIAL

