

PRODUCT INFORMATION

Target	DLL3
Description	Monoclonal Cell Line Derived from 293T Cells, Engineered for Stable Expression of Human DLL3 Using Lentiviral Technology
Host Cells	293T
Uniprot ID	Q9NYJ7
Applications	FACS Data
Growth media	DMEM+10% FBS+1% P.S+Gln+2 ug/mL Puromycin
Package	5E6 Cells/mL
Host Species	Human
Suggested Control	SKU: BME100068
Warranty and Disclaimer	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
Yefei_Storage	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
Synonyms	SCDO1
Background	This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Mutations in this gene cause autosomal recessive spondylocostal dysostosis 1. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]
Usage	For research use only.



Hu_DLL3 293T Cell Line




Figure 1. Flow cytometry analysis of human DLL3 overexpression using Hu_DLL3 293T Cell Line (Cat. No. CEL100039) and Anti-DLL3(Rovalpituzumab biosimilar) mAb (Cat. No. BME100068)

DIMABIO CONFIDENTIAL

