

PRODUCT INFORMATION

Target	GPR75
Description	Monoclonal Cell Line Derived from CHO-S Cells, Engineered for Stable Expression of Human GPR75 Using Lentiviral Technology
Host Cells	CHO-S
Uniprot ID	O95800
Applications	FACS Data
Growth media	DMEM+10% FBS+1% P.S+Gln+2 ug/mL Puromycin
Package	5E6 Cells/mL
Host Species	Human
Suggested Control	SKU: DMC100368
Warranty and Disclaimer	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
Yefei_Storage	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
Synonyms	Probable G-protein coupled receptor 75
Background	GPR75 is a member of the G protein-coupled receptor family. GPRs are cell surface receptors that activate guanine-nucleotide binding proteins upon the binding of a ligand.[supplied by OMIM, Jul 2002]
Usage	For research use only.



Hu_GPR75 N-Strep CHO-S Cell Line




Figure 1. Flow cytometry analysis of Human GPR75 N-Strep overexpression using Hu_GPR75 N-Strep CHO-S Cell Line (Cat. No. CEL100109) and Anti-GPR75 antibody(DMC368); IgG1 Chimeric mAb(Cat. No. DMC100368)

DIMABIO CONFIDENTIAL

