

**PRODUCT INFORMATION**

<b>Target</b>	HVEM
<b>Description</b>	Monoclonal Cell Line Derived from CHO-S Cells, Engineered for Stable Expression of Human HVEM Using Lentiviral Technology
<b>Host Cells</b>	CHO-S
<b>Uniprot ID</b>	Q92956
<b>Applications</b>	FACS Data
<b>Growth media</b>	DMEM+10% FBS+1% P.S+Gln+2 ug/mL Puromycin
<b>Package</b>	5E6 Cells/mL
<b>Host Species</b>	Human
<b>Suggested Control</b>	SKU: DME100131
<b>Warranty and Disclaimer</b>	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
<b>Yefei_Storage</b>	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
<b>Synonyms</b>	ATAR; CD270; HVEA; HVEM; LIGHTR; TR2
<b>Background</b>	This gene encodes a member of the TNF (tumor necrosis factor) receptor superfamily. The encoded protein functions in signal transduction pathways that activate inflammatory and inhibitory T-cell immune response. It binds herpes simplex virus (HSV) viral envelope glycoprotein D (gD); mediating its entry into cells. Alternative splicing results in multiple transcript variants.
<b>Usage</b>	For research use only.



### Hu\_HVEM CHO-S Cell Line




Figure 1. Flow cytometry analysis of human HVEM overexpression using Hu\_HVEM CHO-S Cell Line (Cat. No. CEL100071) and Anti-HVEM antibody(DM131)Rabbit mAb (Cat. No. DME100131)

DIMABIO CONFIDENTIAL

