


**PRODUCT INFORMATION**

<b>Target</b>	SEZ6
<b>Description</b>	Monoclonal Cell Line Derived from K562 Cells, Engineered for Stable Expression of Human SEZ6 Using Lentiviral Technology
<b>Host Cells</b>	K562
<b>Uniprot ID</b>	Q53EL9
<b>Applications</b>	FACS Data
<b>Growth media</b>	RPMI-1640+10% FBS+1% P.S+1% Gln+2 ug/mL Puromycin
<b>Package</b>	5E6 Cells/mL
<b>Host Species</b>	Human
<b>Suggested Control</b>	SKU: DMC101096
<b>Warranty and Disclaimer</b>	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
<b>Storage &amp; Shipping</b>	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
<b>Synonyms</b>	BSRPC
<b>Background</b>	The protein encoded by this gene is thought to contain five cysteine-rich motifs that are similar to sushi domains, as well as two domains similar to the amino terminal half of the CUB (for complement C1r/C1s, Uegf, Bmp1) domain. Mutations in this gene have been associated with febrile seizures. [provided by RefSeq, Jul 2016]
<b>Usage</b>	For research use only.



### Hu\_SEZ6 K562 Cell Line





-  Human IgG
-  Anti-SEZ6 antibody(31A10), IgG1 Chimeric mAb (SKU: DMC101096)

Figure 1. Flow cytometry analysis of human SEZ6 overexpression using Hu\_SEZ6 K562 Cell Line (Cat. No. CEL100037) and Anti-SEZ6 antibody(31A10)IgG1 Chimeric mAb (Cat. No. DMC101096)

DIMABIO CONFIDENTIAL

