

PRODUCT INFORMATION

Target	ADAM8
Synonyms	Cell surface antigen MS2;CD156a
Description	Recombinant human ADAM8 protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	P78325
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	ADAM8(Ile17-Pro655) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 70.6 kDa after removal of the signal peptide. The apparent molecular mass of ADAM8-His is approximately 55-70 kDa and 70-100 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Yefei_Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene may be involved in cell adhesion during neurodegeneration, and it is thought to be a target for allergic respiratory diseases, including asthma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2009]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human ADAM8 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

