

PRODUCT INFORMATION

Tag	C-Flag&Strep Tag
Target	ADGRG2
Synonyms	CBAVDX; EDDM6; GPR64; HE6; TM7LN2
Description	Human ADGRG2-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q8IZP9
Expression Host	HEK293
Protein Families	Druggable Genome, GPCR, Transmembrane
Protein Pathways	N.A.
Molecular Weight	The human full length ADGRG2-Strep protein has a MW of 111.6 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C(Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	A member of the G protein-coupled receptor family described as an epididymis-specific transmembrane protein. The encoded protein may be proteolytically processed as it contains a motif shown to be a protein scission motif in some members of this family.
Usage	Research use only
Conjugate	Unconjugated

