

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C\u00e4wroot\u00fcmirror.dimabio.com\u00e4wp-content\plugins\u00e4woocommerce-print-products\u00e4public\class-woocommerce-print-products-public.php on line 2806 C-flag Tag Tag

AGRB3 Synonyme

Human AGRB3 full length protein-synthetic nanodisc 6~8weeks

Description Delivery Uniprot ID 060242 HEK293 Expression Host

Protein Families GPCR,Transmembrane,Druggable Genome

Protein Pathways Angiogenesis,

Molecular Weight The human full length AGRB3 protein has a MW of 171.5kDa

Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCI, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Store at -20°C to .80°C for 12 omoths in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

This p53-target gene encodes a brain-specific angiogenesis inhibitor, a seven-span transmembrane protein, and is thought to be a member of the secretin receptor family. Brain-specific angiogenesis proteins BAI2 and BAI3 are similar to BAI1 in structure, have similar tissue specificities, and may also play a role in angiogenesis. [provided by RefSeq. Jul 2008]

Research use only

> Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate



