

PRODUCT INFORMATION

Target	C12orf54
Synonyms	HSD-29; HSD-30
Description	Recombinant protein of human chromosome 12 open reading frame 54 (C12orf54)
Delivery	2-3 weeks
Uniprot ID	Q6X4T0
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	14.3 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Recombinant protein was produced with TrueORF clone, [RC206696]. Click on the TrueORF clone link to view cDNA and protein sequences.
Usage	Research use only
Conjugate	Unconjugated

