

PRODUCT INFORMATION

Target	C2orf15
Synonyms	N/A
Description	Recombinant protein of human chromosome 2 open reading frame 15 (C2orf15)
Delivery	2-3 weeks
Uniprot ID	Q8WU43
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	13.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Recombinant protein was produced with TrueORF clone, [RC205732]. Click on the TrueORF clone link to view cDNA and protein sequences.
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

