

PRODUCT INFORMATION

Target	C3orf36
Synonyms	FLJ22173; MGC125760; MGC125761
Description	Purified recombinant protein of Human chromosome 3 open reading frame 36 (C3orf36), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Delivery	2-3 weeks
Uniprot ID	Q3SXR2
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	16.7 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl
Storage & Shipping	Store at -80°C.
Background	A DNA sequence encoding human full-length C3orf36
Usage	Research use only
Conjugate	Unconjugated

DIMA BIO CONFIDENTIAL

