

PRODUCT INFORMATION

Target	C4orf46
Synonyms	RCDG1
Description	Recombinant protein of human chromosome 4 open reading frame 46 (C4orf46)
Delivery	2-3 weeks
Uniprot ID	Q504U0
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	11.7 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	This gene encodes a small, conserved protein of unknown function that is expressed in a variety of tissues. There are pseudogenes for this gene on chromosomes 6, 8, 16, and X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2013]
Usage	Research use only
Conjugate	Unconjugated

