

PRODUCT INFORMATION

Target	CA13
Synonyms	CAXIII
Description	Recombinant protein of human carbonic anhydrase XIII (CA13)
Delivery	2-3 weeks
Uniprot ID	Q8N1Q1
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	29.3 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Carbonic anhydrases (CAs) are a family of zinc metalloenzymes. For background information on the CA family, see MIM 114800.[supplied by OMIM, Mar 2008]
Usage	Research use only
Conjugate	Unconjugated

