

**PRODUCT INFORMATION**

<b>Target</b>	CALML3
<b>Synonyms</b>	CLP
<b>Description</b>	Recombinant protein of human calmodulin-like 3 (CALML3)
<b>Delivery</b>	2-3 weeks
<b>Uniprot ID</b>	P27482
<b>Expression Host</b>	HEK293T
<b>Tag</b>	C-Myc/DDK
<b>Molecular Characterization</b>	N/A
<b>Molecular Weight</b>	16.7 kDa
<b>Purity</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Formulation &amp; Reconstitution</b>	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
<b>Storage &amp; Shipping</b>	Store at -80°C.
<b>Background</b>	May function as a specific light chain of unconventional myosin-10 (MYO10), also enhances MYO10 translation, possibly by acting as a chaperone for the emerging MYO10 heavy chain protein. May compete with calmodulin by binding, with different affinities, to cellular substrates.[UniProtKB/Swiss-Prot Function]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

