

PRODUCT INFORMATION

Target	CCDC112
Synonyms	MBC1
Description	Purified recombinant protein of Human coiled-coil domain containing 112 (CCDC112), transcript variant 2, full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Delivery	1 week
Uniprot ID	Q8NEF3
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	53.4 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150mM NaCl, 10% glycerol, 1% Sarkosyl.
Storage & Shipping	Store at -80°C.
Background	A DNA sequence encoding human full-length CCDC112
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

