

**PRODUCT INFORMATION**

<b>Target</b>	CCDC17
<b>Synonyms</b>	coiled-coil domain containing 17; FLJ17921; FLJ33084; FLJ33084, RP4-697E16.4; OTTHUMP00000009506
<b>Description</b>	Recombinant protein of human coiled-coil domain containing 17 (CCDC17)
<b>Delivery</b>	2-3 weeks
<b>Uniprot ID</b>	Q96LX7
<b>Expression Host</b>	HEK293T
<b>Tag</b>	C-Myc/DDK
<b>Molecular Characterization</b>	N/A
<b>Molecular Weight</b>	48.3
<b>Purity</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Formulation &amp; Reconstitution</b>	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
<b>Storage &amp; Shipping</b>	Store at -80°C.
<b>Background</b>	Recombinant protein was produced with TrueORF clone, [RC206743]. Click on the TrueORF clone link to view cDNA and protein sequences.
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

