

PRODUCT INFORMATION

Target	CCNB1IP1
Synonyms	C14orf18; HEI10
Description	Purified recombinant protein of Human cyclin B1 interacting protein 1, E3 ubiquitin protein ligase (CCNB1IP1), transcript variant 1, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Delivery	2-3 weeks
Uniprot ID	Q9NPC3
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	31.4 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl.
Storage & Shipping	Store at -80°C.
Background	HEI10 is a member of the E3 ubiquitin ligase family and functions in progression of the cell cycle through G(2)/M.[supplied by OMIM, Apr 2004]
Usage	Research use only
Conjugate	Unconjugated

