

PRODUCT INFORMATION

Target CD56

Synonyms NCAM1;CD56;MSK39;NCAM

DescriptionRecombinant Human CD56 Protein with C-

terminal 6×His tag

Delivery In Stock
Uniprot ID P13591
Expression Host HEK293
Tag C-6×His Tag

Molecular Characterization

Molecular Weight

Storage & Shipping

Background

Purity

CD56(Leu20-Gly708) 6×His tag

The protein has a predicted molecular mass of

77.0 kDa after removal of the signal peptide. The apparent molecular mass of CD56-His is

approximately 100-130 kDa due to glycosylation.

The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue

staining.

Lyophilized from sterile PBS, pH 7.4. Normally 5 %
- 8% trehalose is added as protectants before
lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

temperature.

This gene encodes a cell adhesion protein which is a member of the immunoglobulin superfamily. The encoded protein is involved in cell-to-cell interactions as well as cell-matrix interactions during development and differentiation. The

encoded protein plays a role in the development of the nervous system by regulating neurogenesis, neurite outgrowth, and cell migration. This protein is also involved in the expansion of T lymphocytes, B lymphocytes and natural killer (NK) cells which play an important role in immune surveillance. This protein plays a role in signal transduction by interacting with fibroblast growth factor receptors, N-cadherin and other components of the extracellular matrix and by triggering signalling cascades involving FYN-focal adhesion kinase (MAPK), and phosphatidylinosital 3-

protein kinase (MAPK), and phosphatidylinositol 3-kinase (PI3K). One prominent isoform of this gene, cell surface molecule CD56, plays a role in several myeloproliferative disorders such as acute myeloid leukemia and differential expression of this gene is associated with differential disease progression. For example, increased expression of CD56 is correlated with lower survival in acute myeloid leukemia patients whereas increased severity of COVID-19 is correlated with decreased abundance of CD56-expressing NK cells in parisheral blood. Alternative collicing results in

peripheral blood. Alternative splicing results in multiple transcript variants encoding distinct protein isoforms.

Usage Research use only

Conjugate Unconjugated

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China)

/+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com








Figure 1. Human CD56 Protein, His Tag on SDS-PAGE under reducing condition.

Human CD56, His Tagged protein ELISA

0.2 μg of Human CD56, His tagged protein per well




Figure 2. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human CD56 Protein, His Tag (PME100194) can bind Anti-CD56 Neutralizing antibody BME100102 in a linear range of 0.51–41.15 ng/mL.

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

